


## SPECIFICATIONS FID500+Methaniser


Accuracy	Depend of the range used
Drift	1 % over 24 hours
Temperature drift	1 % per degree
Operating temperature	± 20°C without wide variations of temperature
<b>Sampling gas</b>	<b>N<sub>2</sub>, Ar, He, Air, H<sub>2</sub>, O<sub>2</sub></b>
Sample gas connection	1/8" Swagelok SS
Sample flow rate	Approximately 3 to 5 l/h
Sample pressure	Lower than 100 mBar
<b>Combustive gas</b>	<b>Synthetic air</b>
Combustive gas connection	1/8" Swagelok SS
Combustive gas pressure	2 Bar stable
Combustive gas flow rate	300 ml/min
Recommended quality	5.0
<b>Fuel gas</b>	<b>Hydrogen</b>
Fuel gas connection	1/8" Swagelok SS
Fuel gas pressure	1 Bar stable
Fuel gas flow rate	38 ml/min
Recommended quality	6.0
<b>Carrier gas</b>	<b>Argon, Nitrogen, Helium or Hydrogen</b>
Carrier gas connection	1/8" Swagelok SS
Carrier gas pressure	from 4 to 10 Bar
Carrier gas flow rate	2 to 6 l/h
Recommended quality	6.0
Power supply	220 Vac, 50-60 Hz
Power consumption	500 VA
4-20 mA output	Eight configurable outputs depending on the application.
RJ-45 connection	Computerised system maintenance
Output relays (SPST 2 amperes / 250 Vac)	1 Analyser Failure alarm contact 1 Alarm High contact 1 Alarm High High contact

## UMTR Dimensions

Standard UMTR rack mount 2U  
Height > 89 mm | Depth > 365 mm | Width > 483 mm



Standard FID rack mount 5U  
Height > 222 mm | Depth > 545 mm | Width > 483 mm



## FID Dimensions

sales@  
www.orthodyne.be

Line 500



A new generation  
of intelligent detectors

# FID500

+ Methaniser

Analysis of **CH<sub>4</sub>**, **CO**, **CO<sub>2</sub>** and **NMHC**  
in **PPB** and **PPM**



Rue Des Technologies, 23 – B-4432 ALLEUR – BELGIUM  
Phone : +32-4-247.91.06 – Fax : +32-4-263.09.79  
E-Mail : sales@orthodyne.be

**ORTHODYNE**  
GAS CHROMATOGRAPHY

# FID500

+ methaniser

The FID500 is an analytical system that measures CH<sub>4</sub>, CO, CO<sub>2</sub>, NMHC in ppb and in ppm level in Helium, Argon, Oxygen, Nitrogen, Hydrogen or Air

## PRINCIPLE

The flame ionization detector is placed in a temperature regulated chamber.

It is designed to detect traces of hydrocarbons in neutral gases.

The combustion of Hydrogen and Synthetic Air creates a flame in which are burning the organic components contained in the gas to be analysed.

When burning, these components produce ions which are collected by an electrode.

The very weak current obtained in this way is amplified in an electrometer with high gain and directed to a data acquisition system.


A polarization electrode is connected on the level of the nozzle and a collecting electrode with adjustable distance make the best results possible.

Coupled with a methanizer (UMTR unit), it also detects traces of CO and CO<sub>2</sub>.

The methaniser is foreseen to convert, in a catalytic reactor, traces of CO and CO<sub>2</sub> into methane.

This reaction takes place at a temperature of ± 350°C in presence of Hydrogen in excess

The following reactions bring these conversions :




With a properly adjusted temperature and flowrate, the efficiency is almost of 100%.

When supplied with high purity gases, the detection threshold can reach 1 ppb.

## UMTR METHANISER EXPLANATION

- A > Gas inlet
- B > Gas outlet
- C > Oven
- D > Heating element
- E > PT100 sensor
- F > Catalytic reactor




## Type of configuration

- **FID530 (FID530/UMTR)**  
1 Valve /1 column
- **FID540 (FID540/UMTR)**  
1 Valve /2 columns
- **FID550 (FID550/UMTR)**  
2 Valves /1 column
- **FID560 (FID560/UMTR)**  
2 Valves /2 columns
- **FID570 (FID570/UMTR)**  
3 Valves /1 column
- **FID580 (FID580/UMTR)**  
2 Valves /2 columns  
+ external rack

## Applications

- Air separation plants
- Cryogenic truck loading station
- Specialty gas laboratories
- Process control
- Steel industry

## FID DETECTOR EXPLANATION




- A > Sample inlet
- B > Mixture between the sample and the Hydrogen
- C > Hydrogen inlet
- D > Synthetic Air inlet
- E > Nozzle
- F > Flame tip
- G > Collector
- H > Anode & Ignitor
- I > Flame guard

## FEATURES


- < 1 ppb resolution guaranteed. (Limit Detection Level).
- User-friendly software.
- GC technology used for complete separation between each impurity.
- NMHC : Total hydrocarbons from C2 to C5 (Given in CH<sub>4</sub> equivalent).
- Electronic flame-out guard circuit.
- Automatic fuel shut off system.
- Adjustable alarm and oven settings.
- Fast response.
- Possibility of auto-calibration programming – ideal for unmanned plant conditions.
- Easy access to pressure and flow control devices.
- CE marked.

## MEASUREMENT CAPABILITIES

Sample	ARGON	HELIUM	HYDROGEN	NITROGEN	OXYGEN
Type :	< 1 ppb CH <sub>4</sub>	< 1 ppb CH <sub>4</sub>	< 1 ppb CH <sub>4</sub>	< 1 ppb CH <sub>4</sub>	< 1 ppb CH <sub>4</sub>
Orthodyne	< 1 ppb CO	< 1 ppb CO	< 1 ppb CO	< 1 ppb CO	< 1 ppb CO
FID	< 1 ppb CO <sub>2</sub>	< 1 ppb CO <sub>2</sub>	< 1 ppb CO <sub>2</sub>	< 1 ppb CO <sub>2</sub>	< 10 ppb CO <sub>2</sub>
	< 2 ppb NMHC	< 2 ppb NMHC	< 2 ppb NMHC	< 2 ppb NMHC	< 2 ppb NMHC



## System overview



Line 500 - FID500 + Methaniser